

Guidelines

Blood Sample Preparation for Metabolic Phenotyping

Matrix:

Blood samples for quantification of endogenous metabolites using mass spectrometry (metabolomics)

- 1) EDTA plasma* for quantification of endogenous metabolites, e.g. amino acids, biogenic amines, acylcarnitines, phospholipids, hexoses, bile acids, etc.
- 2) EDTA plasma with antioxidant for quantification of eicosanoids
- 3) Serum for quantification of steroid hormones among other uses

Precautions:

- Collect samples in the morning after overnight fasting (before breakfast) or after a fasting period of at least 6 hours prior to sampling.
- When labeling any vials, please ensure the labels are waterproof and resistant to cold storage conditions.
- Please keep all processing procedures and times standardized and use identical blood collection and storage tubes in a single study to ensure comparability.
- When collecting several samples for different analyses, please use the first sample for metabolomics analysis.

Materials:

1) EDTA Plasma

Plasma preparation: S-Monovette 2.7 ml Potassium-EDTA, code red, for plasma separation, with potassium-EDTA, SARSTEDT AG & Co., Nümbrecht, Germany, Art.-No. 05.1167(.001)

Storage vials: Biozym 1.5 ml vial for screw caps, Item no. 710020; Biozym screw cap, transparent, Item no. 710030

Sample volume: Depending on assay requirements

2) EDTA Plasma with antioxidant BHT

Plasma preparation: S-Monovette 2.7 ml Potassium-EDTA, code red, for plasma separation, with potassium-EDTA, SARSTEDT AG & Co., Nümbrecht, Germany, Art.-No. 05.1167(.001)

Storage vials: sampling tubes with BHT (butylated hydroxy toluene), SPI Bio Bertin Pharma, Cat. No. D31007; supplied by Cayman Chemical, Item No. 10950

Sample volume: Depending on assay and sampling tube requirements

3) Serum

Serum preparation: S-Monovette 2.7 ml Z, code white, for serum separation, with additive carrier/clot activator, SARSTEDT AG & Co., Nümbrecht, Germany, Art.-No. 05.1557(.001)

Storage vials: Biozym 1.5 ml vial for screw caps, Item no. 710020; Biozym screw cap, transparent, Item no. 710030

Sample volume: Depending on assay requirements

* Heparin plasma can also be used, Citrate plasma is not recommended

Sample collection, handling and storage (EDTA plasma for amino acids etc.):

- Take blood samples from a peripheral vein directly in tubes for EDTA plasma preparation (see Materials).
- Please ensure the blood sampling tubes are completely filled.
- After collecting blood, shake the tubes gently but thoroughly.
- Do not cool blood before plasma separation has been completed.
- Separate cells and plasma by centrifugation as soon as possible. The time from blood collection to centrifugation should be approximately 40 min. Do not exceed 2 hours. Centrifuge at 20-24 °C for 10 minutes at 2500 x g.
- Transfer plasma into a pre-cooled collection vial (e.g. Falcon) without aspirating blood cells. Use disposable pipette tips; shake plasma thoroughly (Vortex) and place on ice.
- Label the sample storage vials. Cool the sample storage vials and perform the pipetting steps on ice.
- Aliquot plasma into the pre-cooled and labeled storage vials (Biozym, see Materials).
- Freeze plasma aliquots immediately and store at or below -80 °C until shipment. Record the time of collection and the time the samples are placed in the freezer.
- Transport the frozen samples on dry ice according to shipment instructions.

Sample collection, handling and storage (EDTA plasma with BHT for eicosanoids):

- Take blood samples from a peripheral vein directly in tubes for EDTA plasma preparation (see Materials).
- Please ensure the blood sampling tubes are completely filled.
- After collecting blood, shake the tubes gently but thoroughly.
- Do not cool blood before plasma separation has been completed.
- Separate cells and plasma using centrifugation as soon as possible. The time from blood collection to centrifugation should be approximately 40 min. Do not exceed 2 hours. Centrifuge at 20-24 °C for 10 minutes at 2500 x g.
- Transfer plasma into a pre-cooled collection vial (e.g. Falcon) without aspirating blood cells. Use disposable pipette tips; shake plasma thoroughly (Vortex) and place on ice.
- Label the sample storage vials. Cool the sample storage vials and perform the pipetting steps on ice.
- Aliquot plasma into the pre-cooled and labeled sampling tubes with BHT (SPI Bio, Bertin Pharma, see Materials); mix gently but thoroughly.
- Freeze plasma aliquots immediately and store at or below -80 °C until shipment. Record the time of collection and the time the samples are placed in the freezer.
- Transport the frozen samples on dry ice according to shipment instructions.

Sample collection, handling and storage (serum for steroids):

- Take blood samples from a peripheral vein directly in tubes for serum preparation with a clotting activator (see Materials).
- Please ensure the blood sampling tubes are completely filled.
- After collecting blood, shake the tubes gently but thoroughly.
- Store the vial at room temperature (20-24 °C) in upright position to allow coagulation. Clotting is usually completed after 20-30 min. If centrifugation is not performed at the place of sample collection, please use this time for transportation. The time at room temperature until centrifugation should not exceed 40 minutes.
- Centrifuge to separate the serum from the blood clot (15 °C, 10 minutes, 2500 x g).
- Transfer the serum into a pre-cooled collection vial (e.g. Falcon) without aspirating blood cells. Use disposable pipette tips; shake serum thoroughly (Vortex) and place on ice.
- Label the sample storage vials. Cool the sample storage vials and perform the pipetting steps on ice.
- Aliquot serum into the pre-cooled and labeled storage vials (Biozym, see Materials).
- Freeze serum aliquots immediately and store at or below -80 °C until shipment. Record the time of collection and the time the samples are placed in the freezer.
- Transport the frozen samples on dry ice according to shipment instructions.

Sample shipment:

- Please inform the analytical laboratory about the sample shipment 2 to 3 days before the actual shipment.
- Please provide a tracking number.
- Please provide an electronic sample list (use Excel template).
- Package the samples on sufficient dry ice (minimum 10 kg for Europe, 20 kg overseas, thick-walled styrofoam container); the samples should be in labeled boxes protected by a plastic bag.
- The analytical lab will be able to receive samples on working days (8 a.m. to 5 p.m.).

Guidelines

Preparation of Tissue and Feces Samples for Metabolic Phenotyping

Sample collection

- To obtain best results from your metabolomics experiment, prepare all samples in an identical manner as far as possible. If unforeseen deviations occur, make a note of these in your description of the samples. This will help in discussing unexpected results.
- For preparation of tissue or feces samples, prepare 1.5 ml reaction tubes (e.g. 1.5 ml Eppendorf reaction tubes).
- Before transferring the tissue or feces sample into a tube, scale the empty tube and note its weight. This will make it easier to scale the tissue samples.
- Scale the tissue sample. A minimum of 50 mg of tissue is required (in some circumstances, the quantity of tissue has to be adapted to the assay of interest).
- To achieve optimal results, freeze the samples immediately after preparation; otherwise the concentrations of sensitive metabolites may be erroneous.
- Avoid thawing of samples after homogenization and avoid exposure to any organic solvents; this may lead to inaccurate mass spectrometry results. Besides, the samples should not contain detergents.

Homogenization of tissue samples

- Homogenize the tissue sample thoroughly in an appropriate homogenization system (e.g. Precellys Kit or others).
- As extraction buffer we recommend an ethanol/phosphate buffer (3 µl buffer/1 mg tissue).
- Centrifuge the samples for 5 minutes at 5000 x g at 4°C (39.2°F).
- Transfer the supernatant to a tube with the same barcode.
- If the samples are not processed directly after tissue homogenization, all homogenates must be stored immediately at -80°C (-112°F).

Submission of samples

- Submit the samples as soon as possible (preferably on Monday or Tuesday) and send them on dry ice

Sample preparation – Feces samples

- For preparation of feces samples, prepare 1.5 ml reaction tubes (such as 1.5 ml Eppendorf reaction tubes).
- Before transferring the feces sample into a tube, scale the empty tube and note its weight (see tissue sample preparation).
- For sample collection, use a spatula to transfer the feces sample into a 1.5 ml reaction tube and freeze immediately in liquid nitrogen, or transfer the samples to a -80°C (-112°F) freezer.

Homogenization of fecal samples

- For metabolite extractions use an ethanol buffer. The volume of the extraction buffer depends on the metabolites of interest. Mix the samples gently but thoroughly.
- Centrifuge the samples together for 10 minutes at 5000 x g at 4°C (39.2°F).
- Transfer the supernatant to a new labeled 1.5 ml reaction tube and start the analysis.

Sample collection – Urine samples

- Transfer the urine sample in an appropriately labeled 1.5 ml reaction tube. Make sure the tube is locked tight!
- Freeze the sample immediately in liquid nitrogen or on dry ice.
- Avoid freezing-thawing cycles.